

Observasi dengan SEM

By: Yusuf Umardani M.Eng

LPPT-UGM



Outline

- Review perbedaan SEM dan FE-SEM
- Memahami gambar SEM
- **Memilih setting alat yang tepat**
- Bercerita dengan image

SEM vs FE-SEM

Field Emission Gun


- The tip of a tungsten needle is made very sharp (radius < 0.1 μm)
- The electric field at the tip is very strong ($> 10^7$ V/cm) due to the sharp point effect
- Electrons are pulled out from the tip by the strong electric field
- Ultra-high vacuum (better than 10^{-6} Pa) is needed to avoid ion bombardment to the tip from the residual gas.
- Electron probe diameter < 1 nm is possible




Field-emission gun
1st anode to regulate the field strength at the tip by U_1 and hence regulate the emission current.
2nd anode: accelerates the electrons to the final kinetic energy.

LPPT UGM

Source of Electrons



W and LaB₆



Cold- and thermal FEG

Electron Gun Properties

Source	Brightness	Stability(%)	Size	Energy spread	Vacuum
W	3×10^5	~1	50 μm	3.0(eV)	10^{-5} (τ)
LaB ₆	3×10^6	~2	5 μm	1.5	10^{-6}
C-FEG	10^9	~5	5nm	0.3	10^{-10}
T-FEG	10^9	<1	20nm	0.7	10^{-9}

Brightness – beam current density per unit solid angle

LPPT UGM

Membaca image SEM

Contoh: Struktur keramik

Informasi pada image:

- Jenis detector
- Tegangan yg dipakai
- Working distance
- Spot size
- Perbesaran
- Scaler (mistar bar).

LPPT UGM

Kualitas gambar berbeda jika level tegangan berbeda

Effect of acceleration voltage

High resolution can be obtained.

High

↑

Acceleration voltage

↓

Low

Resolution lowers.

Image of surface structure is sharp. Damage is less. Charge-up is reduced. Edge effect is reduced.

Image of surface structure is dull. Edge effect is large. Charge-up easily occurs. Damage is more significant.

Diffusion of incident electrons (by Duncumb and Shields)

Low acceleration voltage

High acceleration voltage

Lower atomic No.

Higher atomic No.

LPPT UGM

Effect of Accelerating Voltage

Specimen: Toner

When high accelerating voltage is used as at (a), it is hard to obtain the contrast of the specimen **surface structure**. Besides, the specimen surface is easily charged up. The surface microstructures are easily seen at (b).

<http://lppt.uqm.ac.id>


Effect of Accelerating Voltage

Specimen: Evaporated Au particles

The image sharpness and resolution are better at the higher accelerating voltage, 25 kV.


<http://lppt.uqm.ac.id>

Efek Probe Current




LPPT UGM

Efek Working Distance (WD) pada image




LPPT UGM

Diamete aperture terhadap image



LPPT UGM

Influence of Working Distance and Aperture Size




<http://ppt.ugm.ac.id>

Influence of Working Distance and Aperture Size



(d) OL aperture diameter: 200µm WD: 38mm





(e) OL aperture diameter: 100µm WD: 38mm

Specimen Electric bulb coil
5 kV x540
The smaller the OL aperture diameter and the longer the WD, the greater the depth of field.

<http://lppt.ugm.ac.id>

Pengamatan Sample non conductive



LPPT UGM

Tip dan trik sample

- Menentukan jenis sample: untuk memahami preparasi yg cocok
- Menentukan titik yang akan diamati: permukaan atau cross section
- Konduktif atau non konduktif
- Mengatur setting alat
- Pengaturan khusus utk EDS dan mapping
- Analisis menyeluruh

LPPT UGM

Topography and morphology

- High depth of focus






Image: Christian Kjølseth, UiO


Image: Camilla Kongshaug, UiO

BSE vs SE




Images: Greg Meeker, USGS

Bercerita dengan gambar



LPPT UGM

Contoh sample lempung



LPPT UGM